

FORM MD-5

[See sub-rule (4) of rule 20 and sub-rule (6) of rule 20]

Licence to Manufacture for Sale or for Distribution of Class A or Class B medical device

Licence Number: MFG/MD/2024/000906

- 1. M/s HIPPOCRAT BIOMEDICAL LLP, 1402 NEW MHADA COMPLEX CLUSTER 4 L15 MIRA ROAD KASHIMIRA POLICE ST VASAI THANE 401107 MAHARASHTRAMIRA ROAD, Thane, Maharashtra (India) 401107 Telephone No.: 9323982285 FAX: 9136652285 has been licenced to manufacture for sale or for distribution the below listed medical device(s) at the premises situated at M/s HIPPOCRAT BIOMEDICAL LLP, GALA NO C/10, Unit No DV01/1347, GANESH INDUSTRIAL ESTATE, DEVDAL, SURVEY NUMBER 43,HISSA NO 1/2/B, NEAR SAGAR HOTEL, BHIWANDI ROAD, KAMAN, VASAI, Palghar, Maharashtra (India) 401208 Telephone No.: 9323982285 FAX: 9323982285
- 2. Details of medical device(s) [Annexed]
- 3. This licence is subject to the provisions of the Medical Devices Rules, 2017 and conditions prescribed therein.

ANNEXURE

S.No.	Details Of Device(s)
1	Generic Name:Hydrogen Peroxide and Peracetic Acid Based Cold Sterilant Model No.:HBML01 - Hydrogen peroxide(20%) Peracetic Acid(4.5%) Intended Use:Hydrogen Peroxide and Peracetic Acid-based cold sterilant is intended for cleaning and disinfection of dialysis consumables such as dialyzer and dialysis tubing.
	Class of medical device:Class B Material of construction:Hydrogen Peroxide, Peracetic Acid
	Dimension(if any):NA
	Shelflife:12 Months Sterile or Non sterile:Non-Sterilized
	Brand Name(if registered under the Trade Marks Act, 1999):HIPPOCLEAN
2	Generic Name:Sodium Hypochlorite Model No.:HBML03 - Sodium Hypochlorite (5%)
	Intended Use:Sodium Hypochlorite intended to disinfect hospital equipment Class of medical device:Class B
	Material of construction:Sodium Hypochlorite
	Dimension(if any):NA
	Shelflife:12 Months
	Sterile or Non sterile:Non-Sterilized
	Brand Name(if registered under the Trade Marks Act, 1999):HYPOCLEAN
3	Generic Name:Citric Acid based Hot Disinfectant
	Model No.:HBML04 - Citric Acid 21%,HBML05 - Citric Acid 50% Intended Use:Citric Acid-based Hot Disinfectant intended for the disinfection of
	dialysis machines
	Class of medical device:Class B
	Material of construction:Citric acid
	Dimension(if any):NA
	Shelflife:12 Months
	Sterile or Non sterile:Non-Sterilized
	Brand Name(if registered under the Trade Marks Act, 1999):HIPPOSTERI

